

Human CLIC4 (His tag) recombinant protein



Catalog Number: 518398

General Information

Gene Name Synonym

Chloride Intracellular Channel Protein 4,
Intracellular Chloride Ion Channel Protein p64H1,
CLIC4

Protein Construction

A DNA sequence encoding the Human CLIC4
Met1-Lys253 is expressed with a 6His tag at the
N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM
TrisHCl, 100mM NaCl, 1mM DTT, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CLIC4 (His tag) recombinant protein