

# CLOCK antibody

Catalog Number: 109398

## Product name

CLOCK antibody

IP: 1:200-1:2000

IF: 1:20-1:200

## Specificity

Human; other species not tested.

## Antibody description

CLOCK Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HeLa cells, HEK-293 cells, PC-3 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 150 kDa

## Preparation

This antibody was obtained by immunization of CL<sup>o</sup>CK recombinant protein (Accession Number: XM\_047416433). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IP, IF

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109398(CLOCK Antibody) at dilution of 1:1000



IP Result of anti-CLOCK (IP:Catalog No:109398, 3ug; Detection:Catalog No:109398 1:500) with HEK-293 cells lysate 1000ug.



Immunofluorescent analysis of HeLa cells using Catalog No:109398(CLOCK Antibody) at dilution of

# CLOCK antibody



Catalog Number: 109398

---

1:50 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)