

# Anti-IL-3 antibody



Catalog Number: 100798

## Product name

Polyclonal

Anti-IL-3 antibody

## Ig Type

## Immunogen

Rabbit IgG

[Rat IL-3 \(His Tag\) recombinant protein](#)

## Applications

ELISA, IHC-P

## Specificity

Rat IL-3

## Dilutions

ELISA: 0.1-0.2 µg/mL

## Antibody description

Rabbit polyclonal to IL-3

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat IL-3. The detection limit for Rat IL-3 is approximately 0.00975 ng/well.

## Preparation

Produced in rabbits immunized with purified, recombinant Rat IL-3 (rR IL-3; NP\_113701.1; Met1-Cys169). IL-3 specific IgG was purified by Rat IL-3 affinity chromatography.

IHC-P: 0.1-2 µg/mL

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Validations


## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality



IL-3 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of rat IL3 in rat spleen with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).