

# BNIP3L Polyclonal Antibody



Catalog Number: 164798

## Product name

BNIP3L Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to BNIP3L

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human BNIP3L (NP\_004322.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations




Western blot - BNIP3L Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using BNIP3L antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.



Immunohistochemistry - BNIP3L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using BNIP3L antibody at dilution



of 1:100 (40x lens).

# BNIP3L Polyclonal Antibody




Catalog Number: 164798

---

## Immunohistochemistry - BNIP3L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using BNIP3L antibody at dilution of 1:100



(40x lens).

## Immunohistochemistry - BNIP3L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using BNIP3L antibody at dilution of 1:100 (40x lens).