

# RPAP2 antibody

Catalog Number: 114798

## Product name

RPAP2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RPAP2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, MCF7 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 70kd

## Preparation

This antibody was obtained by immunization of RPAP2 recombinant protein (Accession Number: NM\_024813). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114798(RPAP2 antibody) at dilution of 1:500



IP Result of anti-RPAP2 (IP:Catalog No:114798, 4ug; Detection:Catalog No:114798 1:1000) with HEK-293 cells lysate 2800ug.