

# ADORA2A antibody



Catalog Number: 107798

## Product name

ADORA2A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ADORA2A Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue.

## Preparation

This antibody was obtained by immunization of ADORA2A recombinant protein (Accession Number: NM\_001278498). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:107798(ADORA2A Antibody) at dilution of 1:100 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:107798(ADORA2A Antibody) at dilution of 1:100 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).