

H3N1 HA (His Tag) recombinant protein



Catalog Number: 504001

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/swine/Korea/PZ72-1/2006(H3N1)) hemagglutinin (ACS71642.1) (Met1-Arg345), termed as HA1, was expressed with a polyhistidine tag at the C-terminus.

Organism

H3N1

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The recombinant HA1 subunit of the Influenza A virus hemagglutinin

(A/swine/Korea/PZ72-1/2006(H3N1)) consists 340 amino acids and predicts a molecular mass of 38 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H3N1 (A/swine/Korea/PZ72-1/2006)] Hemagglutinin / HA1 Protein (His Tag)