

ZNF512 antibody



Catalog Number: 117001

Product name

ZNF512 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF512 Rabbit Polyclonal antibody. Positive WB detected in mouse small intestine tissue, A549 cells, fetal human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 66-70 kDa

Preparation

This antibody was obtained by immunization of ZNF512 recombinant protein (Accession Number: NM_001271288). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:117001(ZNF512 Antibody) at dilution of 1:600



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:117001(ZNF512 Antibody) at dilution of 1:300