

Anti-TAF3 Rabbit antibody



Catalog Number: 169001

Product name

Anti-TAF3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TAF3 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human TAF3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

1:1000

1:50-1:200

Validations



Western blot detection of TAF3 in Jurkat cell lysates using TAF3 Rabbit mAb(1:1000 diluted).Predicted band size:104kDa.Observed band size:140kDa.