

# BACE1 antibody

Catalog Number: 108401

## Product name

BACE1 antibody

## Specificity

Human, Mouse, Zebrafish; other species not tested.

## Antibody description

BACE1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse brain tissue, HeLa cells, human brain tissue, mouse pancreas tissue, zebrafish tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 70-75 kDa

## Preparation

This antibody was obtained by immunization of BACE1 recombinant protein (Accession Number: NM\_012104). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP, IHC

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108401(BACE1 antibody) at dilution of 1:500




IP Result of anti-BACE1 (IP:Catalog No:108401, 3ug; Detection:Catalog No:108401 1:500) with mouse brain tissue lysate 7500ug.

# BACE1 antibody

Catalog Number: 108401

---



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:108401(BACE1 Antibody) at dilution of 1:50 (under 40x lens)