

Human UBE2B/RAD6B (His tag) recombinant protein



Catalog Number: 519401

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 B, RAD6 Homolog B, HR6B, hHR6B, Ubiquitin Carrier Protein B, Ubiquitin-Conjugating Enzyme E2-17 kDa, Ubiquitin-Protein Ligase B, UBE2B, RAD6B

Protein Construction

A DNA sequence encoding the Human UBE2B/RAD6B Met1-Ser152 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM PB,150mM NaCl,pH7.2.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2B/RAD6B (His tag) recombinant protein