

RPE antibody

Catalog Number: 114801

Product name

RPE antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RPE Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Positive IF detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of RPE recombinant protein (Accession Number: NM_199229). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114801(RPE antibody) at dilution of 1:200



Immunofluorescent analysis of RAW 264.7 cells using Catalog No:114801(RPE Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)