

VSTM1 antibody

Catalog Number: 116801

Product name

VSTM1 antibody

Specificity

Human; other species not tested.

Antibody description

VSTM1 Rabbit Polyclonal antibody. Positive FC detected in HL-60 cells. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of VSTM1 recombinant protein (Accession Number: XM_017026666). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116801(VSTM1 Antibody) at dilution of 1:300



1X10⁶ HL-60 cells were stained with 0.2ug VSTM1 antibody (Catalog No:116801, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.