

DNAJC1 antibody

Catalog Number: 110101

Product name

DNAJC1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNAJC1 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in Raji cells, C6 cells, NIH/3T3 cells. Positive IP detected in Raji cells. Observed molecular weight by Western-blot: 64kd

Validations




Preparation

This antibody was obtained by immunization of DNAJC1 recombinant protein (Accession Number: BC030955). Purification method: Antigen affinity purified.

Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:110101(DNAJC1 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110101(DNAJC1 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:2000



IP Result of anti-DNAJC1 (IP:Catalog No:110101,

DNAJC1 antibody



Catalog Number: 110101

3ug; Detection:Catalog No:110101 1:300) with Raji cells lysate 2000ug.