

GMPPA antibody



Catalog Number: 111101

Product name

GMPPA antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GMPPA Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IF detected in Hela. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of GMPPA recombinant protein (Accession Number: XM_005246486). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111101(GMPPA antibody) at dilution of 1:1000



Immunofluorescent analysis of Hela cells, using GMPPA antibody Catalog No:111101 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).