

Anti-MSC antibody



Catalog Number: 176101

Product name

Anti-MSC antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to MSC

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human musculin. aa range:131-180

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1.; Western Blot analysis of various cells using Musculin Polyclonal Antibody diluted at 1: 2000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).