

TMEM161A antibody



Catalog Number: 116101

Product name

TMEM161A antibody

Specificity

Human; other species not tested.

Antibody description

TMEM161A Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells, HT-1080 cells, PC-3 cells. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of TMEM161A recombinant protein (Accession Number: NM_017814). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116101(TMEM161A Antibody) at dilution of



1:1000

MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116101(TMEM161A Antibody) at dilution of 1:1000