

Sabatolimab (Synonyms: ANTI-TIM3 CHECKPOINT INHIBITOR MBG453)



Catalog Number: 177101

Product name

Sabatolimab (Synonyms: ANTI-TIM3 CHECKPOINT INHIBITOR MBG453)

Specificity

Human TIM-3 / HAVCR2 / CD366

Antibody description

Anti-TIM-3 / HAVCR2 / CD366 Antibody (Sabatolimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4SP

Applications


ELISA

Validations



SDS-PAGE

Anti-TIM-3 / HAVCR2 / CD366 Antibody (sabatolimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TIM-3 / HAVCR2 / CD366 Antibody (sabatolimab) is 98.49%, determined by SEC-HPLC.