

# Anti-Cytochrome C (7C10) Mouse antibody



Catalog Number: 168101

## Product name

Mouse IgG1

Anti-Cytochrome C (7C10) Mouse antibody

## Applications

## Specificity

WB, ICC/IF, IHC-P

Human, Mouse, Rat, Chicken

## Dilutions

## Antibody description

WB: 1/1000-5000

Cytochrome C (7C10) Mouse monoclonal antibody

IHC: 1/500-1000

## Preparation

IF: 1/200

Antigen: Recombinant Protein of CYCS

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage


Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

## Validations



Western blot detection of Cytochrome c in T47D, 3T3, COS7, C6 and HeLa cell lysates using Cytochrome c mouse mAb (1:10000 diluted). Predicted band size: 14 kDa. Observed band size: 14 kDa.