

## General Information

### Gene Name Synonym

Immunoglobulin heavy chain gamma polypeptide

### Protein Construction

A DNA sequence encoding the Mouse IgG2a Fc region (P01863) (Glu98-Lys330) was expressed and purified.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 92 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Glu 98

### Molecular Mass

The recombinant mouse IgG2a-Fc consists of 233

amino acids and has a predicted molecular mass of 26 kDa. As a result of glycosylation, the apparent molecular mass of IgG2a-Fc is approximately 34 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse IgG2a-Fc Protein SDS-PAGE