

EXOC1 antibody



Catalog Number: 110501

Product name

EXOC1 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

IP: 1:500-1:5000

[Human EXOC1 Recombinant protein \(GST tag\)](#)

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

EXOC1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in rat brain tissue, human brain tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 102 kDa




Preparation

This antibody was obtained by immunization of EXOC1 recombinant protein (Accession Number: XM_017008409). Purification method: Antigen affinity purified.

rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110501(EXoc1 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110501(EXoc1 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB, IP

Dilutions

EXOC1 antibody



Catalog Number: 110501



IP Result of anti-EXOC1 (IP:Catalog No:110501, 3ug; Detection:Catalog No:110501 1:1000) with mouse brain tissue lysate 3600ug.