

# TYROBP Polyclonal Antibody



Catalog Number: 164501

## Product name

TYROBP Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to TYROBP

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human TYROBP (NP\_003323.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - TYROBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using TYROBP antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - TYROBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using TYROBP antibody at dilution of 1:100 (40x lens).