

# Anti-KCTD7 antibody



Catalog Number: 176501

## Product name

Anti-KCTD7 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to KCTD7

## Preparation

This antigen of this antibody was synthesized peptide derived from human kctd7. at aa range: 181-230

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1:; Western Blot analysis of 1,mouse-brain 2,mouse-spleen cells using primary antibody diluted at 1:1000(4°C overnight). Secondary antibody: Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25°C, 1 hour)