

# RBM4 antibody



Catalog Number: 107501

## Product name

RBM4 antibody

## Specificity

Human; other species not tested.

## Antibody description

RBM4 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, fetal human brain tissue. Positive IHC detected in human cervical cancer tissue, human stomach tissue. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of RBM4 recombinant protein (Accession Number: BC021120). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107501(RBM4 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:107501(RBM4 Antibody) at dilution of 1:50 (under 40x lens)