

KLRK1 Polyclonal Antibody



Catalog Number: 160901

Product name

KLRK1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to KLRK1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human KLRK1 (NP_031386.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - KLRK1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using KLRK1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.