

Muskelin antibody

Catalog Number: 112901

Product name

Muskelin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Muskelin Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, Jurkat cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 70-75 kDa

Preparation

This antibody was obtained by immunization of Muskelin recombinant protein (Accession Number: NM_013255). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:112901(MKLN1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112901(MKLN1 antibody) at dilution of 1:50 (under 40x lens)