

Anti-JAK2 Rabbit antibody



Catalog Number: 168901

Product name

Anti-JAK2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

JAK2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human JAK2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

1:1000

1:20

Validations



Western blot detection of JAK2 in K562,3T3 cell lysates using JAK2 Rabbit mAb(1:1000 diluted).Predicted band size:131kDa.Observed band size:131kDa.