

# PDC Polyclonal Antibody



Catalog Number: 162201

---

## Product name

PDC Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to PDC

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human PDC (NP\_002588.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Validations



Immunohistochemistry - PDC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat retina using PDC antibody at dilution of 1:200 (40x lens).