

# SERPINB12 antibody



Catalog Number: 115201

## Product name

SERPINB12 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

SERPINB12 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 46-48kd

## Preparation

This antibody was obtained by immunization of SERPINB12 recombinant protein (Accession Number: NM\_001307928). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115201(SERPINB12 antibody) at dilution of 1:500