

Anti-CD19 Mouse mAb (FITC)



Catalog Number: 168201

Product name

Anti-CD19 Mouse mAb (FITC)

Specificity

Human

Antibody description

CD19 Mouse mAb (FITC)

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Mouse IgG1

Applications

ICC/IF

Dilutions

ICC: 1/100

Validations

Confocal immunofluorescent analysis of Ramos (positive cell, left) and Jurkat (negative cell, right) using anti-CD19-FITC mouse mAb (dilution 1:100).