

# Anti-CD19 Mouse mAb (FITC)



Catalog Number: 168201

---

**Product name**

Mouse IgG1

Anti-CD19 Mouse mAb (FITC)

**Applications****Specificity**

ICC/IF

Human

**Dilutions****Antibody description**

ICC: 1/100

CD19 Mouse mAb (FITC)

**Validations****Formulation**

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

**Storage**


Store at 4°C short term.Avoid freeze / thaw cycle.

**Clonality**


Monoclonal

**Ig Type**

Ramos



Jurkat



Confocal immunofluorescent analysis of Ramos (positive cell, left) and Jurkat (negative cell, right) using anti-CD19-FITC mouse mAb (dilution 1:100).