

# RNASE3 Polyclonal Antibody



Catalog Number: 163601

## Product name

RNASE3 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to RNASE3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 28-160 of human RNASE3 (NP\_002926.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:100

## Validations



Immunohistochemistry - RNASE3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using RNASE3 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - RNASE3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using RNASE3 antibody at dilution of 1:100 (40x lens).