

USP24 antibody

Catalog Number: 116601

Product name

USP24 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

USP24 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 211 kDa

Preparation

This antibody was obtained by immunization of USP24 recombinant protein (Accession Number: BC029660). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116601(USP24 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using USP24 antibody Catalog No:116601 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).