

Anti-Pyruvate Dehydrogenase E2 Rabbit antibody



Catalog Number: 167601

Product name

Anti-Pyruvate Dehydrogenase E2 Rabbit antibody

Applications

WB, ICC/IF, FC, IP

Specificity

Human, Mouse, Rat

Dilutions

WB: 1/1000

Antibody description

Pyruvate Dehydrogenase E2 Rabbit monoclonal antibody

ICC/IF: 1/20

FC: 1/20

Preparation

Antigen: A synthetic peptide of human Pyruvate Dehydrogenase E2

IP: 1/20

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Western blot detection of Pyruvate Dehydrogenase E2 in CHO-K1,C6,Rat Brain,Ramos cell lysates using Pyruvate Dehydrogenase E2 Rabbit mAb(1:1000 diluted).Predicted band size:69kDa.Observed band size:69kDa.

