

SRP9 antibody



Catalog Number: 107601

Product name

SRP9 antibody

Specificity

Human; other species not tested.

Antibody description

SRP9 Mouse Monoclonal antibody. Positive WB detected in Hela cells. Observed molecular weight by Western-blot: 10kd

Preparation

This antibody was obtained by immunization of SRP9 recombinant protein (Accession Number: NM_003133). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107601(SRP9 antibody) at dilution of 1:500