

Anti-Activin A Receptor Type IB Rabbit antibody



Catalog Number: 168601

Product name	Rabbit IgG
Anti-Activin A Receptor Type IB Rabbit antibody	Applications
Specificity	WB, IHC-P, ICC/IF
Human, Mouse, Rat	Dilutions
Antibody description	WB: 1/1000-1/5000
Activin A Receptor Type IB Rabbit monoclonal antibody	IHC: 1/50-1/200
Preparation	ICC/IF: 1/20-1/50

Antigen: A synthetic peptide of human Activin A Receptor Type IB

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Validations



Western blot detection of Activin A Receptor Type IB in CHO-K1,Rat Brain lysates using Activin A Receptor Type IB Rabbit mAb(1:500 diluted).Predicted band size:57kDa.Observed band size:57kDa.