

Human LGALS1 (His tag) recombinant protein



Catalog Number: 519601

General Information

Gene Name Synonym

Galectin-1, Gal-1, 14 kDa Laminin-Binding Protein, HLBP14, 14 kDa Lectin, Beta-Galactoside-Binding Lectin L-14-I, Galaptin, HBL, HPL, Lactose-Binding Lectin 1, Lectin Galactoside-Binding Soluble 1, Putative MAPK-Activating Protein PM12, S-Lac Lectin 1, LGALS1

Protein Construction

A DNA sequence encoding the Human LGALS1 Ala2-Asp135 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2 \mu\text{m}$ filtered solution of 10mM PB, 200mM NaCl, 2mM DTT, pH 7.0.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\mu\text{g/ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human LGALS1 (His tag) recombinant protein