

# RPS6KA4 Polyclonal Antibody



Catalog Number: 160301

## Product name

RPS6KA4 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to RPS6KA4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human RPS6KA4 (NP\_003933.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:200 - 1:1000

IHC 1:20 - 1:100


IF1:20 - 1:100

## Validations



Western blot - RPS6KA4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using RPS6KA4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry - RPS6KA4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using RPS6KA4 antibody at dilution of 1:200 (40x lens).