

HESX1 antibody

Catalog Number: 111301

Product name

HESX1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HESX1 Rabbit Polyclonal antibody. Positive IHC detected in human brain.

Preparation

This antibody was obtained by immunization of HESX1 recombinant protein (Accession Number: XM_005265526). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human brain using Catalog No:111301(HESX1 antibody) at dilution of 1:50 (under 10x lens)