

# LPCAT2 antibody



Catalog Number: 112301

## Product name

LPCAT2 antibody

## Immunogen

[Human LPCAT2 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

LPCAT2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 53kd

## Preparation

This antibody was obtained by immunization of LPCAT2 recombinant protein (Accession Number: NM\_017839). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112301(LPCAT2 antibody) at dilution of 1:400