

Mouse anti RFP-Tag mAb



Catalog Number: 166301

Product name

Mouse anti RFP-Tag mAb

Antibody description

Monoclonal antibody to Mouse anti RFP-Tag

Preparation

Antigen: Recombinant protein of human RFP-Tag

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB, IHC, IP


Dilutions

WB 1:2000 - 1:5000

IHC 1:50 - 1:200

IP 1:50 - 1:100

Validations



Western blot - Mouse anti RFP-Tag mAb

Western blot analysis of recombinant RFP protein with RFP-tag antibody at dilution of 1:10000. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.