

TRIM17 antibody



Catalog Number: 116301

Product name

TRIM17 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRIM17 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, COLO 320 cells. Positive IF detected in COLO 320 cells. Observed molecular weight by Western-blot: 54kd

Preparation

This antibody was obtained by immunization of TRIM17 recombinant protein (Accession Number: XM_011544210). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116301(TRIM17 antibody) at dilution of 1:300



Immunofluorescent analysis of COLO 320 cells using Catalog No:116301(TRIM17 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)