

FOXL2-Specific antibody



Catalog Number: 110701

Product name

FOXL2-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FOXL2-Specific Rabbit Polyclonal antibody.
Positive WB detected in mouse ovary tissue.
Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_023067).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:110701(FOXL2-Specific Antibody) at dilution of 1:300