

# HSPB7 antibody



Catalog Number: 111701

## Product name

HSPB7 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

HSPB7 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HeLa cells, HepG2 cells, human brain tissue, human heart tissue. Positive IP detected in mouse skeletal muscle tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 23 kDa

## Validations




## Preparation

This antibody was obtained by immunization of HSPB7 recombinant protein (Accession Number: NM\_014424). Purification method: Antigen affinity purified.

mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111701(HSPB7 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111701(HSPB7 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:



# HSPB7 antibody

Catalog Number: 111701

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111701(HSPB7 antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-HSPB7 (IP:Catalog No:111701, 4ug; Detection:Catalog No:111701 1:500) with mouse skeletal muscle tissue lysate 3600ug.