

RBM6 Monoclonal Antibody



Catalog Number: 164701

Product name

RBM6 Monoclonal Antibody

Specificity

Human, Mouse

Antibody description

Monoclonal antibody to RBM6

Preparation

Antigen: A synthetic peptide of human RBM6

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - RBM6 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using RBM6 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.