

# STBD1 antibody



Catalog Number: 115701

## Product name

STBD1 antibody

## Immunogen

[Human STBD1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

STBD1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, A549 cells, HeLa cells, human heart tissue. Positive IP detected in A549 cells. Observed molecular weight by Western-blot: 35-43kd

## Preparation

This antibody was obtained by immunization of STBD1 recombinant protein (Accession Number: NM\_003943). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:115701(STBD1 antibody) at dilution of 1:1000



IP Result of anti-STBD1 (IP:Catalog No:115701, 4ug; Detection:Catalog No:115701 1:1000) with A549 cells lysate 1600ug.