

V5-tag antibody

Catalog Number: 116701

Product name

V5-tag antibody

Specificity

Human; other species not tested.

Antibody description

V5-tag Rabbit Polyclonal antibody. Positive IP detected in Recombinant protein protein. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of Peptide. Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IP: 1:5000-1:50000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:117342(V5-tag antibody) at dilution of 1:32000



IP Result of anti-V5-tag (IP:Catalog No:117342, 3ug; Detection:Catalog No:117342 1:20000) with Recombinant protein protein lysate 3000ug.