

Anti-Phospholipase C gamma 1/PLC-gamma-1 Rabbit antibody



Catalog Number: 167701

Product name

Rabbit IgG

Anti-Phospholipase C gamma 1/PLC-gamma-1 Rabbit antibody

Applications

WB, ICC/IF, FC, IP

Specificity

Human, Mouse, Rat

Dilutions

WB: 1/2000

Antibody description

Phospholipase C gamma 1/PLC-gamma-1 Rabbit monoclonal antibody

ICC/IF: 1/50

FC: 1/20

Preparation

Antigen: A synthetic peptide of human Phospholipase C gamma 1/PLC-gamma-1

IP: 1/20

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Phospholipase C gamma 1/PLC-gamma-1 in U2OS,C6 using Phospholipase C gamma 1/PLC-gamma-1 Rabbit mAb(1:1000 diluted)