

DR3 antibody

Catalog Number: 110011

Product name

DR3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DR3 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse spleen tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 17-20 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003790). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110011(TNFRSF25 antibody) at dilution of 1:600