

KIF2B antibody

Catalog Number: 112011

Product name

KIF2B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KIF2B Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse brain tissue, mouse testis tissue. Observed molecular weight by Western-blot: 80 kDa

Preparation

This antibody was obtained by immunization of KIF2B recombinant protein (Accession Number: NM_032559). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112011(KIF2B antibody) at dilution of 1:1000



IP Result of anti-KIF2B (IP:Catalog No:112011, 4ug; Detection:Catalog No:112011 1:1000) with mouse testis tissue lysate 4000ug.